



# HbA1C

Awareness



## HbA1C

by HPLC Method (BIORAD D-10)



• A0

• HbA1c

• Total area 1M-4M

• A1c result

### Uses of HbA1C

- ❖ Gold Standard Methodology.
- ❖ Fasting sample not required.
- ❖ Persistent Hyperglycemia is better measured by HbA1C
- ❖ Monitoring the efficiency of Dietary Control and therapy



### Advantages of HPLC Method

- ❖ Excellent specificity and reproducibility.
- ❖ Separates Hb species based on charge and detect more common HB Variable.
- ❖ Not influenced by low Hb / Heterophile Antibody
- ❖ Affected in haematological disorder such as polcythemia and incorrect sample collection

Kindly ask for HbA1C Graph when you visit Lab / Hospital

Check for articles on google by typing search word - HbA1C by HPLC method

Still got doubts, talk to our friendly pathologist for more information about tests done in our Lab

## BHAGAVAN PATHOLOGY LABORATORY

Main Center : # 1116, 5th Cross, 1503, Srirampet, Vinoba Road | ☎ 0821-4256488, 2424488

Branch : # 665, 17th Main, Kamakshi Hospital Road, Mysuru-09 | ☎ 0821 - 4269488

🌐 www.bhagavanpathlab.com | ✉ bhagavanpathlab@hotmail.com | reports@bhagavanpathlab.com

Whatsapp Helpline No. : 📞 8762424488